 PAUL HARTMANN AG · 89522 Heidenheim, Germany

ES — Laboratorios HARTMANN S.A. · 08302 Mataró
GR — PAUL HARTMANN Hellas A.E. · 16674 Glyfada /Athina
PT — PAUL HARTMANN LDA · 2685-378 Prior Velho
ZA — HARTMANN South Africa · 2194 Johannesburg

www.hartmann.info

helps healing.

925 082/0 (0511)

Thermoval® duo scan

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

1. Preliminary remarks

Please read these instructions for use carefully before first use as correct temperature measurement depends only on the appropriate use of the device. These instructions describe the individual steps of temperature measurement using Thermoval duo scan and contain important and helpful hints for the reliable determination of the body temperature. Be sure to keep these instructions for use for future reference.

Table of contents

	Page
1. Preliminary remarks	43
2. Introduction	44
3. Benefits of the thermometer	44
4. General information on body temperature	44
5. Safety precautions	45
6. Initial operation of the device	45
7. Measurement of the body temperature	45
7.1. Ear measurement	46
7.2. Forehead measurement	47
8. Cleaning and maintenance of the device	49
9. Power supply and disposal notes	49
9.1. Batteries and disposal	49
9.2. Battery replacement	49
10. Error messages	50
11. Calibration check	50
12. Warranty conditions	50
13. Signs and symbols	50
14. Technical data	51
15. Legal requirements and guidelines	51
16. Contact information for customer queries	52

2. Introduction

Dear Customer, we are pleased that you have decided to purchase a clinical thermometer from HARTMANN. Thermoval duo scan is a quality product for measuring the human body temperature in the ear and on the forehead. It is ideally suited for measurements on children as from 6 months of age, but the thermometer can be used also on adults. Given correct application, the device ensures a fast and precise measurement of the body temperature in a very comfortable manner. We wish you all the best for your health.

3. Benefits of the thermometer

Two measuring methods: Ear and forehead measurement

Thermoval duo scan measures the temperature in the ear and on the forehead. Depending on your individual preference, you can choose between these two modes. The device recognizes automatically whether an ear measurement or forehead measurement is carried out.

High measuring accuracy ensured by the latest measuring sensor technology

The latest generation of the measuring sensors guarantees a very high measuring accuracy in ear and forehead measurements.

Quick measurement ensured by infrared technology

The measuring sensor picks up the infrared radiation emitted by the eardrum and/or the surface of the forehead and displays the measuring result within a second (three seconds in case of forehead measurement).

Easy to use

The device is very easy and safe to use having been provided with two buttons only. The O/I button is used to switch the device on (or off), and the SCAN button starts the measurement. In addition, the illustrations on the storage box explain how to carry out the measurements correctly.

Robust and reliable use ensured by high quality manufacturing

This device has been developed on the basis of the strict HARTMANN quality guidelines. Geared to the customer requirements, it fulfils the highest standards in stability.

4. General information on body temperature

The human body adjusts the body temperature to a setpoint value – it varies up to 1 °C throughout the day. In addition, the temperature in the body (core temperature) and the surface temperature of the skin show different temperatures. For this reason, there is no “normal” body temperature – it always depends on the measurement site. Here as well, the body temperature is influenced by the ambient temperature, age, stress, duration of sleep, hormones and physical activity. Whilst glass thermometers and digital thermometers measure the temperature of the human body directly, thermometers using ear and forehead measurement determine the core temperature through the infrared radiation of the body. This measurement may slightly deviate from the temperature of rectal, oral or axillary measurements using digital thermometers.

Measurement table (in °C):

Designation	Measurement site		
	Ear / forehead	Anus	Mouth / armpit
Subnormal temperature	< 35.7	< 36.2	< 35.8
Normal temperature	35.8 – 36.9	36.3 – 37.4	35.9 – 37.0
Elevated temperature	37.0 – 37.5	37.5 – 38.0	37.1 – 37.5
Slight fever	37.6 – 38.0	38.1 – 38.5	37.6 – 38.0
Moderate fever	38.1 – 38.5	38.6 – 39.0	38.1 – 38.5
High fever	38.6 – 39.4	39.1 – 39.9	38.6 – 39.5
Very high fever	39.5 – 42.0	40.0 – 42.5	39.6 – 42.0

5. Safety precautions

- This device may be used only for the measurement of human body temperature in the ear and on the forehead. The manufacturer shall not be held liable for any damage resulting from improper handling.

- Do not leave the device unattended near toddlers or persons who cannot operate it themselves.
- Do not expose the device to hard knocks or vibrations. Do not use the device, if you notice any kind of damage.
- Prior to every measurement, please check whether the measuring sensor is clean and unscratched.
- Please observe the storage and operating conditions mentioned in Chapter “Technical data”. Protect the device against dirt and dust, extreme temperatures, solar radiation and fluids.
- For cleaning, please follow the instructions mentioned in Chapter “Cleaning and maintenance of the device”.
- The device must not be altered, dismantled or repaired by the user.
- The use does not replace medical treatment in any way.
- If you inform your doctor about the temperature measured, please also mention that this temperature has been determined in the ear or on the forehead.

6. Initial operation of the device

Batteries are included in the scope of delivery and have been placed in the device already. However, a contact strip interrupts the contact of the batteries to the device and thus guarantees the full battery capacity prior to the first use of the thermometer. Please carefully pull out the contact strip from the closed lid of the battery compartment. Subsequently the Thermoval duo scan is ready for operation.

7. Measurement of the body temperature

Thermoval duo scan measures the temperature on the forehead as well as in the ear. Both types of

measurement can be used for the measurement of your own body temperature (self-measurement) as well as for the measurement of the body temperature of another person (foreign measurement).

1. Before starting the measurement, the device has to be switched on. For this reason, press the ON/OFF button (O/I). All symbols will be shown on the display for two seconds.
2. The device displays the temperature measured last for three seconds. At the same time you will see a symbol in the display indicating whether an ear measurement  or a forehead measurement  has been carried out.
3. Subsequently the temperature display will disappear and the thermometer symbol  will be faded in. Now the device is ready for measurement:

Important information to prevent inaccurate measured values:


- Please observe that the clinical thermometer and the person to be measured are in a room at room temperature for at least 30 minutes prior to the measurement.
- In order to obtain the highest of accuracy in successive measurements, please wait for 30 seconds between two measurements.
- After sleeping, it is recommended to wait for a few minutes before measuring the temperature.
- Do not measure temperature directly after showering, swimming, et cetera, when the ears or the forehead are wet.
- Avoid eating, drinking or exercising prior to the measurement.

- Do not measure the temperature during or directly after breast feeding a baby.
- Remove the measuring device from the measuring area only after the end signal has been sounded and the LED lamp is no longer lit.
- In the following situations, we recommend to carry out three temperature measurements and to consider the highest value as measuring result:
 - Children below the age of three with a weak immune system (in particular, if the presence or absence of fever is significant).
 - In case of users, who are not familiar with the device, until they receive identical measured values.
 - In case of suspiciously low measured values.


7.1 Ear measurement

The measurement can be carried out in the right-hand or left-hand ear. The thermometer measures the infrared rays emitted by the eardrum and the ear canal. These rays are picked up by the sensor and converted into temperature values. The best values are obtained when the infrared rays are picked up from the middle ear (eardrum). Rays which are emitted from surrounding tissue produce lower measured values.


If the cap for forehead measurement has been placed, please remove it by simply pulling it. For better grip, both sides of the cap have been provided with rough surfaces.


Stretch the ear canal by pulling the middle of the outer ear to the back and up at the same time. In case of children below the age of one year, pull the outer ear to the back only. Carefully insert the thermometer tip into the ear canal until you feel a slight resistance.

Press the SCAN button to start measurement. At the beginning of measurement a brief tone sounds.

During the duration of the measurement, the LED over the SCAN button is lit. The ear measurement takes one second. When the measurement is finished, a long tone is sounded and the measuring result is displayed.

Please clean the thermometer tip directly after the measurement using the cleaning swabs enclosed.

The thermometer symbol  will appear three seconds after the measurement, and the device is ready for the next measurement.

Important information for ear measurement:

- If earwax accumulates on the measuring sensor, this may lead to inaccuracies of the temperature measurement or may even cause an infection between various users. For this reason, it is imperative that the measuring sensor is cleaned for every measurement. For notes on cleaning, please read the instructions in Chapter "Cleaning and maintenance of the device".
- After cleaning the measuring sensor with alcohol-impregnated swabs, please wait for 5 minutes until carrying out the next measurement to permit the thermometer to reach the necessary operating temperature.
- The temperature values may vary from ear to ear. For this reason, please always use the same ear


for measuring the temperature.

- No disposable protective caps are required for the measuring sensor. Please do not use disposable protective caps as they may distort the measuring result.
- If you have rested on an ear, the temperature in this ear may be elevated shortly after.
- The measurement must not be carried out in an ear which is suffering from inflammatory diseases as well as after ear injuries (such as eardrum injury) or during the healing phase after a surgical procedure. If medication has been applied to an ear, this ear must not be used for measurement either. In these cases, please consult your attending doctor.


7.2 Forehead measurement

The thermometer measures the infrared radiation emitted by the skin in the forehead and temple area and the surrounding tissue. The radiation is picked up by the sensor and converted into temperature values. The most precise values are obtained when you scan the area over the eyebrows towards the temple.


If the cap for forehead measurement has not been placed, please place it by simply pressing it on again.

Position the surface of the forehead adapter about 1 cm over the middle of the eyebrows.

Press the SCAN button to start measurement. At the beginning of measurement a brief tone sounds.

Move the thermometer along the forehead towards the temple. During the duration of the measurement, the LED lamp over the SCAN button is lit. At the end of measurement, a long tone is sounded and the measuring result is displayed. The forehead measurement takes three seconds.

48

The thermometer symbol  will appear three seconds after the measurement, and the device is ready for the next measurement.

Important information for forehead measurement:

- We recommend to measure the temperature at the same position every time as otherwise the displayed values may vary.
- Remove hair, perspiration, cosmetics or dirt from the forehead to ensure an optimum measurement accuracy.

Switching over from ear to forehead measurement (or vice versa) is effected automatically by means of a contact sensor in the thermometer tip. The thermometer recognizes automatically which mode is operational.

8. Cleaning and maintenance of the device

The measuring sensor is the most important and most sensitive part of the device. In order to ensure most precise measurements, it always has to be clean and intact. Clean the measuring tip and the sensor exclusively with the Soft-Zellin alcohol swabs supplied along with the device or with a soft moist cloth or with 70 percent alcohol (isopropanol). Residues from earwax can be removed carefully using cotton buds. Please do not use any aggressive cleaning agents or solvents as they may damage the device and could make the display intransparent. Please

English **GB**

make sure that no fluid can penetrate the interior of the thermometer. To protect it from external influences keep the thermometer in the storage box.

9. Power supply and disposal notes

9.1 Batteries and disposal

- The two high-quality batteries included with the device guarantee at least 1,000 measurements. Use only high-quality batteries (see specification in Chapter "Technical data").
- If you are using less efficient batteries, we can no longer guarantee 1,000 measurements.
- Never mix old and new batteries or batteries made by different manufacturers.
- Immediately remove exhausted batteries.
- If you do not intend to use the device for a longer period of time, you should remove the batteries to prevent possible leakage.
- Please help protect the environment: Do not dispose of batteries in household waste! Use designated collection points or municipal collection / recycling centres to dispose of hazardous household waste.

This symbol on products and/or accompanying documents means that spent electronic products must not be mixed with regular household waste.

9.2 Battery replacement

This device has been provided with two alkaline batteries of LR03 (AAA) type. Please replace spent

49

batteries by two new LR03 batteries as soon as the following battery symbol is indicated in the display:  .

Open the battery cover in the direction of the arrow on the rear side of the device by applying gentle pressure on the notch. Now insert the two “AAA” size batteries so that the positive (+) and negative (–) polarities match the “+” and “–” signs shown on the marking in the battery compartment. If the polarities do not match, the device will not work and the batteries may leak! Replace the lid of the battery compartment, applying gentle pressure until the hook snaps in place. The last measured value recorded in the memory will not be lost when changing the batteries.

10. Error messages

Thermoval duo scan is a clinically validated, high-quality product. Nonetheless error messages may occur, such as in the case of excessively high or low ambient temperature, if the measuring temperature is outside the range for human body temperatures, if the batteries are empty or in very rare cases, if a system error occurs. The description of errors is provided at the end of these instructions for use along with the explanation of the symbols on the display.

11. Calibration check

We recommend a calibration check at intervals of two years in the case of professionally used devices, for example in pharmacies, medical practices or hospitals. You should also observe the national regulations determined by the legislator. Calibration checks can be carried out either by the competent authorities or authorized maintenance providers against compensation.

12. Warranty conditions

We give a **3-year warranty** on this high-quality clinical ear and forehead thermometer from the date of purchase and in accordance with the following conditions. Claims must be made during the warranty period. The date of purchase may be documented by the appropriately completed and stamped warranty certificate or proof of purchase.

Within the warranty period, HARTMANN shall replace or repair any faulty device components free of charge which were caused by material or manufacturing errors. This does not extend the warranty period.

This warranty is not applicable to damage caused by improper use or unauthorized interference. Parts that are subject to wear and tear (batteries, forehead adapter, storage box, packaging et cetera) are excluded from the warranty. Claims for compensation are limited to the value of the goods; compensation for subsequent damages is expressly excluded.

In warranty cases please send the device with cap, storage box and with the fully completed and stamped warranty certificate direct, or via your dealer to the Customer Services department for your country.

AU – PAUL HARTMANN Pty. Ltd.
Level 6, 5 Ryder Boulevard
Rhodes, NSW 2138, Australia

ZA – HARTMANN South Africa
2194 Johannesburg

13. Signs and symbols

 Disposal note for electronic devices

 Disposal note

 Observe the operating instructions

 Please note

 Protection from electric shock (type BF)

 PAUL HARTMANN AG
89522 Heidenheim, Germany

 0 1 2 3

14. Technical data

Type:	Infrared thermometer Thermoval duo scan
Measurement range:	Ear mode: 32.0 °C – 42.2 °C Forehead mode: 34.0 °C – 42.2 °C
Measuring accuracy:	Laboratory: ± 0.2 °C at 35.5 °C – 42.0 °C ± 0.3 °C at 32.0 / 34.0 °C – 35.4 °C and 42.1 °C – 42.2 °C in accordance with European Standard EN 12470-5
Display:	Liquid crystal display, 4 digits plus special symbols
Smallest display unit:	0.1 °C
Acoustic signal:	1 short signal tone when measuring process starts 1 long signal tone when measuring process ends
Optical signal:	Blue LED lamp is lit during the measuring process
Memory:	Last measured value
Background illumination:	The display is lit in blue when the device has been switched on and after measuring processes
Operating conditions:	Ambient temperature: 16 °C – 35 °C Relative humidity: 30 % – 85 %

Storage / transport conditions:	Ambient temperature: -25 °C – 55 °C Relative humidity: 30 % – 85 % Humidity for storage: 15 % – 85 %
Automatic switch-off:	About 1 minute after the last measurement
Power supply:	2 x 1.5 V Alkaline-manganese Mignon (AAA/LR03)
Battery capacity:	At least 1,000 measurements
Size:	150 mm (l) x 38 mm (w) x 40 mm (h)
Weight:	About 100 g (with battery)
Serial number:	See battery compartment

15. Legal requirements and guidelines

Thermoval duo scan complies with the requirements of the EC directive 93/42/EEC on medical devices (Medical Device Directive MDD) and bears the CE mark.

European Standard DIN EN 12470-5: Performance of infrared ear thermometers
ISO 80601-2-56

Portable and mobile high-frequency and communication devices can impair the functional capability of electronic medical devices.

16. Contact information for customer queries

AU – PAUL HARTMANN Pty. Ltd.
Level 6, 5 Ryder Boulevard
Rhodes, NSW 2138
Australia

ZA – HARTMANN South Africa
2194 Johannesburg

Date of revision of the text: 2011-03